

Anti-HLA-A, B, C antibody [W6/32] {EF700} (STJA0016417)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	EF700
Purification	Affinity purified
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID [3106](#)
[3105](#)
[3107](#)
[HLA-B](#)
[HLA-A](#)
[HLAB_HUMAN](#)
[HLAA_HUMAN](#)
[HLAC_HUMAN](#)

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.