

**Anti-HLA-A, B, C antibody [W6/32] {EV450} (STJA0016418)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC                 |
| Host / Source | Mouse              |
| Reactivity    | Human              |

**PRODUCT PROPERTIES**

|                     |                                                                                                  |
|---------------------|--------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                       |
| Concentration       | 5 Mu L/Test                                                                                      |
| Conjugation         | EV450                                                                                            |
| Purification        | Affinity purified                                                                                |
| Formulation         | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.      |
| Isotype             | IgG2ak                                                                                           |
| Storage Instruction | Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

**TARGET INFORMATION**

Gene ID [3106](#)  
[3105](#)  
[3107](#)  
[HLA-B](#)  
[HLA-A](#)  
[HLAB\\_HUMAN](#)  
[HLAA\\_HUMAN](#)  
[HLAC\\_HUMAN](#)

**ADDITIONAL INFORMATION**

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.