

Anti-CD11b antibody [ICRF44] {EV500} (STJA0016426)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	FC
	Host / Source	Mouse
	Reactivity	Human
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Conjugation	EV500
	Purification	Affinity purified
	Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
	Isotype	IgG1k
	Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.
TARGET INFORMATION		
	Gene ID	3684
	Gene Symbol	ITGAM
	UniProt ID	ITAM_HUMAN
ADDITIONAL INFORMATION		
	Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.