

Anti-CD3 antibody [G4.18] {EF700} (STJA0016476)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	EF700
Purification	Affinity purified
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG3k
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	300678 25710
Gene Symbol	Cd3g Cd3d
UniProt ID	CD3G_RAT CD3D_RAT

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---