

**Anti-CD158b/j antibody [DX27] {PerCP/Cyanine 5.5} (STJA0016694)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	5 Mu L/Tests
Conjugation	PerCP/Cyanine 5.5
Purification	Affinity purified
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

**TARGET INFORMATION**

Gene ID	<a href="#">3804</a> <a href="#">3803</a>
Gene Symbol	<a href="#">KIR2DL3</a> <a href="#">KIR2DL2</a>
UniProt ID	<a href="#">K12L3_HUMAN</a> <a href="#">K12L2_HUMAN</a>

**ADDITIONAL INFORMATION**

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---