

## Anti-CD23 antibody [EBVCS2] {PerCP/Cyanine 5.5} (STJA0016699)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5 Mu L/Tests
Conjugation	PerCP/Cyanine 5.5
Purification	Affinity purified
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

### TARGET INFORMATION

Gene ID	<a href="#">2208</a>
Gene Symbol	<a href="#">FCER2</a>
UniProt ID	<a href="#">FCER2_HUMAN</a>

### ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------