

Anti-CD192/CCR2 antibody [K036C2] {APC} (STJA0016708)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | FC |
| Host / Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Conjugation | APC |
| Purification | Affinity purified |
| Formulation | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant. |
| Isotype | IgG2ak |
| Storage Instruction | Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

| | |
|-------------|----------------------------|
| Gene ID | 729230 |
| Gene Symbol | CCR2 |
| UniProt ID | CCR2_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|---|