

Anti-CD3/CD4/CD8a/CD45 Cocktail antibody [OKT-3;RPA-T4;OKT-8;HI30] {FITC/PE/PE/Cyanine 7/PerCP} (STJA0016757)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	FITC/PE/PE/Cyanine 7/PerCP
Purification	Affinity purified
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---