

Anti-Recombinant-beta Actin antibody [SM8244] {HRP} (STJA0018244)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	HRP
Purification	Affinity purification
Dilution Range	WB (H, M, R) 1:1000-1:3000
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Storage Instruction	Store at-20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11461
Gene Symbol	Actb
UniProt ID	ACTB_MOUSE

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Sample: Protein treated by RIPA Lysis Buffer.
 Blocking buffer: 3% Nonfat dry milk in TBST, RT,
 1 h. Antibody: HRP conjugated anti-beta actin
 antibody 1:3000, 4C overnight..