

## Anti-CD2 antibody [B-E2] {R-PE} (STJA0018250)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	R-PE
Purification	Protein A purification
Dilution Range	Please carry out application specific titration to determine the most appropriate dilution ratio.
Formulation	PBS with 0.2% BSA, and 0.09% Sodium Azide
Isotype	IgG2b
Storage Instruction	Shipped on gel ice packs. Store at 2-8C on arrival.

### TARGET INFORMATION

Gene ID	<a href="#">914</a>
Gene Symbol	<a href="#">CD2</a>
UniProt ID	<a href="#">CD2_HUMAN</a>

### ADDITIONAL INFORMATION

Note	FOR RESEARCH USE ONLY (RUO).
------	------------------------------