

Anti-CD4 antibody [Edu-2] {CyQ} (STJA0018252)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	CyQ
Purification	Protein A purification
Dilution Range	Please carry out application specific titration to determine the most appropriate dilution ratio.
Formulation	PBS with 0.2% BSA, and 0.09% Sodium Azide
Isotype	IgG2a
Storage Instruction	Shipped on gel ice packs. Store at 2-8C on arrival.

TARGET INFORMATION

Gene ID	920
Gene Symbol	CD4
UniProt ID	CD4_HUMAN

ADDITIONAL INFORMATION

Note	FOR RESEARCH USE ONLY (RUO).
------	------------------------------