

Anti-CD15 antibody [BRA-4F1] {FITC} (STJA0018275)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	FITC
Purification	Protein A purification
Dilution Range	Please carry out application specific titration to determine the most appropriate dilution ratio.
Formulation	PBS with 0.2% BSA, and 0.09% Sodium Azide
Isotype	IgM
Storage Instruction	Shipped on gel ice packs. Store at 2-8C on arrival.

ADDITIONAL INFORMATION

Note FOR RESEARCH USE ONLY (RUO).

