

## Anti-CD16 antibody [B-E16] (STJA0018281)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A purification
Dilution Range	Please carry out application specific titration to determine the most appropriate dilution ratio.
Formulation	PBS with 0.2% BSA, and 0.09% Sodium Azide
Isotype	IgG2a
Storage Instruction	Shipped on gel ice packs. Store at 2-8C on arrival.

### TARGET INFORMATION

Gene ID	<a href="#">2214</a> <a href="#">2215</a>
Gene Symbol	<a href="#">FCGR3A</a> <a href="#">FCGR3B</a>
UniProt ID	<a href="#">FCG3A_HUMAN</a> <a href="#">FCG3B_HUMAN</a>

### ADDITIONAL INFORMATION

Note	FOR RESEARCH USE ONLY (RUO).
------	------------------------------