

Anti-CD23 antibody [B-G6] (STJA0018300)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A purification
Dilution Range	Please carry out application specific titration to determine the most appropriate dilution ratio.
Formulation	PBS with 0.2% BSA, and 0.09% Sodium Azide
Isotype	IgG1
Storage Instruction	Shipped on gel ice packs. Store at 2-8C on arrival.

TARGET INFORMATION

Gene ID	2208
Gene Symbol	FCER2
UniProt ID	FCER2_HUMAN

ADDITIONAL INFORMATION

Note	FOR RESEARCH USE ONLY (RUO).
------	------------------------------

