

Anti-CD45RA antibody [MB1] {R-PE} (STJA0018327)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	R-PE
Purification	Protein A purification
Dilution Range	Please carry out application specific titration to determine the most appropriate dilution ratio.
Formulation	PBS with 0.2% BSA, and 0.09% Sodium Azide
Isotype	IgG1
Storage Instruction	Shipped on gel ice packs. Store at 2-8C on arrival.

TARGET INFORMATION

Gene ID	5788
Gene Symbol	PTPRC
UniProt ID	PTPRC_HUMAN

ADDITIONAL INFORMATION

Note	FOR RESEARCH USE ONLY (RUO).
------	------------------------------