

## Anti-CD45RB antibody [MT4] (STJA0018328)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC                 |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Clonality           | Monoclonal  |
| Conjugation         | Unconjugated  |
| Purification        | Protein A purification  |
| Dilution Range      | Please carry out application specific titration to determine the most appropriate dilution ratio. |
| Formulation         | PBS with 0.2% BSA, and 0.09% Sodium Azide   |
| Isotype             | IgG1  |
| Storage Instruction | Shipped on gel ice packs. Store at 2-8C on arrival.   |

### TARGET INFORMATION

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">5788</a>        |
| Gene Symbol | <a href="#">PTPRC</a>       |
| UniProt ID  | <a href="#">PTPRC_HUMAN</a> |

### ADDITIONAL INFORMATION

Note FOR RESEARCH USE ONLY (RUO).

