

Anti-CD103 antibody [B-ly7] {FITC} (STJA0018347)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	FITC
Purification	Protein A purification
Dilution Range	Please carry out application specific titration to determine the most appropriate dilution ratio.
Formulation	PBS with 0.2% BSA, and 0.09% Sodium Azide
Isotype	IgG1
Storage Instruction	Shipped on gel ice packs. Store at 2-8C on arrival.

TARGET INFORMATION

Gene ID	3683
Gene Symbol	ITGAL
UniProt ID	ITAL_HUMAN

ADDITIONAL INFORMATION

Note	FOR RESEARCH USE ONLY (RUO).
------	------------------------------

