

Anti-MPO antibody [266.6K2] {APC} (STJA0018370)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | FC |
| Host / Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Monoclonal |
| Conjugation | APC |
| Purification | Protein A purification |
| Dilution Range | Please carry out application specific titration to determine the most appropriate dilution ratio. |
| Formulation | PBS with 0.2% BSA, and 0.09% Sodium Azide |
| Isotype | IgG1 |
| Storage Instruction | Shipped on gel ice packs. Store at 2-8C on arrival. |

TARGET INFORMATION

| | |
|-------------|----------------------------|
| Gene ID | 4353 |
| Gene Symbol | MPO |
| UniProt ID | PERM_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|------------------------------|
| Note | FOR RESEARCH USE ONLY (RUO). |
|------|------------------------------|

