

Anti-CMV pp65 antibody [C11] (STJA0018377)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Dilution Range	Please carry out application specific titration to determine the most appropriate dilution ratio.
Formulation	Cell culture supernatant with 0.09% Sodium azide, and Fetal Calf Serum
Isotype	IgG1
Storage Instruction	Shipped on gel ice packs. Store at 2-8C on arrival.

ADDITIONAL INFORMATION

Note	FOR RESEARCH USE ONLY (RUO).
------	------------------------------