

Mouse IgG1 Isotype Control [MCG1] {Unconjugated} (STJA0018382)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A purification
Dilution Range	Please carry out application specific titration to determine the most appropriate dilution ratio.
Formulation	PBS with 0.2% BSA, and 0.09% Sodium Azide
Isotype	IgG1
Storage Instruction	Shipped on gel ice packs. Store at 2-8C on arrival.

ADDITIONAL INFORMATION

Note	FOR RESEARCH USE ONLY (RUO).
------	------------------------------