

## Anti-Keratin 8 antibody [C51] (STJA0018397)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Dilution Range	Please carry out application specific titration to determine the most appropriate dilution ratio.
Formulation	Cell culture supernatant with 0.09% Sodium azide, and Fetal Calf Serum
Storage Instruction	Shipped on gel ice packs. Store at 2-8C on arrival.

### TARGET INFORMATION

Gene ID	<a href="#">3856</a>
Gene Symbol	<a href="#">KRT8</a>
UniProt ID	<a href="#">K2C8_HUMAN</a>

### ADDITIONAL INFORMATION

Note	FOR RESEARCH USE ONLY (RUO).
------	------------------------------