

## Anti-IL6 antibody (Val30-Met212) [10G10] (STJA0018412)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | ELISA/IHC/WB       |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Clonality           | Monoclonal  |
| Concentration       | 3.65 mg/mL  |
| Conjugation         | Unconjugated  |
| Purification        | Protein A/G purified from cell culture supernatant  |
| Formulation         | 0.01M PBS, pH 7.4.  |
| Isotype             | IgG1  |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months |

### TARGET INFORMATION

|                  |                           |
|------------------|---------------------------|
| Gene ID          | <a href="#">3569</a>      |
| Gene Symbol      | <a href="#">IL6</a>       |
| UniProt ID       | <a href="#">IL6_HUMAN</a> |
| Immunogen Region | Val30-Met212              |

### ADDITIONAL INFORMATION

|      |  |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS |
|------|--|