

Anti-CDKN2A/p16INK4a antibody (Met1-Asp156) [2E5] (STJA0018415)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IHC/ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Protein G purified from mice ascities
Formulation	0.01M PBS
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

TARGET INFORMATION

Gene ID	1029
Gene Symbol	CDKN2A
UniProt ID	CDN2A_HUMAN
Immunogen Region	Met1-Asp156

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--