

Anti-Abeta42/Beta-APP42 antibody [21F12] (STJA0018426)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IHC/WB
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1.47 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

TARGET INFORMATION

Gene ID	351
Gene Symbol	APP
UniProt ID	A4_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--