

Anti-Recombinant-6*His Tag antibody (HHHHHH) [3D5] (STJA0018431)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF/IP/WB
Host / Source	Rabbit

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.83 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

TARGET INFORMATION

Immunogen Region	HHHHHH
------------------	--------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--