

**Anti-HRSV F/Fusion glycoprotein F0 antibody [Am14] (STJA0018433)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/IHC/WB
Host / Source	Human
Reactivity	HRSV-A

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	1.67 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--