

Anti-HLA-A2-peptide Complex antibody (SLLMWITQV) [3M4E5] (STJA0018478)**GENERAL INFORMATION**

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | ELISA/FC/SPR |
| Host / Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Monoclonal |
| Concentration | 2.19 mg/mL |
| Conjugation | Unconjugated |
| Purification | Protein A/G purified from cell culture supernatant |
| Formulation | 0.01M PBS, pH 7.4. |
| Isotype | IgG1 |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months |

TARGET INFORMATION

| | |
|------------------|-----------|
| Immunogen Region | SLLMWITQV |
|------------------|-----------|

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS |
|------|--|