

Anti-IgA antibody {FITC} (STJA0018527)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IF/WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	FITC
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Isotype	IgG1
Storage Instruction	Store at 4°C 12 months

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--