

**Anti-HRSV-A/BRSV F/Fusion glycoprotein F0-strain A2 & Rb94 antibody [B4-1] (STJA0018560)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/IF/RIA/WB
Host / Source	Human
Reactivity	HRSV-A/BRSV

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

**ADDITIONAL INFORMATION**

Note	<b>STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS</b>
------	---