

Anti-HPV16 E6/Protein E6 antibody [Nb9] (STJA0018595)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/Neut/WB
Host / Source	Mouse
Reactivity	Human papillomavirus type 16

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	2.5 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	VHH-mFC
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--

