

Anti-HCV NS1/gp68/gp70/Envelope glycoprotein E2 antibody [HC84-26] (STJA0018606)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/Neut
Host / Source	Human
Reactivity	HCV

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--