

**Anti-HCV NS1/gp68/gp70/Envelope glycoprotein E2-strain H77 antibody [mAb#12]
(STJA0018629)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/Neut
Host / Source	Human
Reactivity	HCV/Recombinant Hepatitis C virus H77C/JFH-1

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--