

**Anti-CHIKV p130/Structural polyprotein antibody [DC2112] (STJA0018657)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	BIF/ELISA/FC/IP
Host / Source	Human
Reactivity	CHIKV

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	0.84 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--