

Anti-RABV G/Glycoprotein antibody [RVC20] (STJA0018678)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC/Neut
Host / Source	Mouse
Reactivity	Rabies lyssavirus

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	2.27 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--