

## Anti-RABV G/Glycoprotein-strain CVS-11 antibody [Fab094] (STJA0018681)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IF/IP
Host / Source	Mouse
Reactivity	Rabies virus

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1.97 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--