

Anti-DNA/RNA G-quadruplex Structures antibody [BG4] (STJA0018701)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/EMSA/IF/ChIP-seq/WB
Reactivity	General

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1.25 mg/mL
Conjugation	Unconjugated
Purification	Purified by Immobilized Metal Affinity Chromatography
Formulation	0.01M PBS, pH 7.4.
Isotype	scFv-FLAG-His
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--------------------------------------------------------------------------------------------------------------------