

## Anti-SKB1/PRMT5 IP antibody [SAA0228] (STJA0018743)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IF/IP
Host / Source	Mouse
Reactivity	Human/Mouse

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

### TARGET INFORMATION

Gene ID	<a href="#">10419</a>
Gene Symbol	<a href="#">PRMT5</a>
UniProt ID	<a href="#">ANM5_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--