

## Anti-HMGB1 antibody [DPH11] (STJA0019032)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/FuncS/IF/WB
Host / Source	Rat
Reactivity	Human/Mouse/Rat/Bovine/Dog/Horse/Pig

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

### TARGET INFORMATION

Gene ID	<a href="#">15289</a>
Gene Symbol	<a href="#">Hmgb1</a>
UniProt ID	<a href="#">HMGB1_MOUSE</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--