

## Anti-TP53/p53 antibody [PAb242] (STJA0019035)

### GENERAL INFORMATION

|               |                        |
|---------------|------------------------|
| Product Type  | Primary antibodies     |
| Applications  | ELISA/IF/IHC/IP/RIA/WB |
| Host / Source | Rat                    |
| Reactivity    | Mouse                  |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Clonality           | Monoclonal  |
| Concentration       | 1 mg/mL   |
| Conjugation         | Unconjugated  |
| Purification        | Protein A/G purified from cell culture supernatant  |
| Formulation         | 0.01M PBS, pH 7.4.  |
| Isotype             | IgG2a   |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months |

### TARGET INFORMATION

|             |                           |
|-------------|---------------------------|
| Gene ID     | <a href="#">22059</a>     |
| Gene Symbol | <a href="#">Tp53</a>      |
| UniProt ID  | <a href="#">P53_MOUSE</a> |

### ADDITIONAL INFORMATION

|      |  |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS |
|------|--|