

Anti-CXCL10/IP-10 antibody [54] (STJA0019044)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC
Host / Source	Rat
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.75 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

TARGET INFORMATION

Gene ID	18186
Gene Symbol	Nrp1
UniProt ID	NRP1_MOUSE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--