

## Anti-Flag Tag antibody (DYKDDDDK) [M2] (STJA0019061)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ChIP/ELISA/EIA/FC/ICC/IF/IP/WB
Host / Source	Mouse

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	2.6 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

### TARGET INFORMATION

Immunogen Region	DYKDDDDK
------------------	----------

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--------------------------------------------------------------------------------------------------------------------