

Anti-HA Tag antibody (YPYDVPDYA) [12CA5] (STJA0019101)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF/IP/WB
Host / Source	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	2.7 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG2b
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

TARGET INFORMATION

Immunogen Region	YPYDVPDYA
------------------	-----------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	---