

Anti-V5 Tag antibody (IPNPLLGLD) [PK1] (STJA0019110)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IF/IHC-F/IP/RIA/WB
Host / Source	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG2a
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

TARGET INFORMATION

Immunogen Region	IPNPLLGLD
------------------	-----------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--