

Anti-EE Tag antibody (EYMPME) [3D5/EE_48K] (STJA0019114)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/SPR/WB
Host / Source	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1.9 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

TARGET INFORMATION

Immunogen Region	EYMPME
------------------	--------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--