

Anti-ALFA Tag antibody (SRLEEELRRRLTE) [NbALFA] (STJA0019115)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC/IF/IP/SPR/WB
Host / Source	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1.42 mg/mL
Conjugation	Unconjugated
Purification	Protein A or G purified
Formulation	0.01M PBS
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

TARGET INFORMATION

Immunogen Region	SRLEEELRRRLTE
------------------	---------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--