

Anti-Phospho-Gamma-H2AX-S139 antibody [G2] (STJA0019154)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | IF/IHC/IP/WB |
| Host / Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Monoclonal |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Purification | Protein A/G purified from cell culture supernatant |
| Formulation | 0.01M PBS, pH 7.4. |
| Isotype | IgG2a |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months |

TARGET INFORMATION

| | |
|-------------|----------------------------|
| Gene ID | 3014 |
| Gene Symbol | H2AX |
| UniProt ID | H2AX_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS |
|------|--|