

**Anti-Phospho-POLR2A-CTD antibody (YSPTSPSYSPTSPS) [1stS5] (STJA0019159)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IHC/IP/WB
Host / Source	Mouse
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG2a
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

**TARGET INFORMATION**

Gene ID	<a href="#">5430</a>
Gene Symbol	<a href="#">POLR2A</a>
UniProt ID	<a href="#">RPB1_HUMAN</a>
Immunogen Region	YSPTSPSYSPTSPS

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--