

## Anti-KIRREL3/NEPH2 antibody [N321C/49R] (STJA0019224)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF/IHC/WB
Host / Source	Mouse
Reactivity	Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG2a
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

### TARGET INFORMATION

Gene ID	<a href="#">67703</a>
Gene Symbol	<a href="#">Kirrel3</a>
UniProt ID	<a href="#">KIRR3_MOUSE</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--