

## Anti-K48-linked Polyubiquitin antibody [SAA0329] (STJA0019256)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IP/WB
Host / Source	Mouse
Reactivity	Human/Mouse/Rat/Rabbit

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1.08 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Isotype	IgG2a
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

### TARGET INFORMATION

Gene ID	<a href="#">7316</a>
Gene Symbol	<a href="#">UBC</a>
UniProt ID	<a href="#">UBC_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--