

Anti-MAPT/Tau/PHF-tau antibody [KW1] (STJA0019261)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Blocking/ELISA
Host / Source	Human
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Formulation	0.01M PBS, pH 7.4.
Isotype	VHH-Fc
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

TARGET INFORMATION

Gene ID	4137
Gene Symbol	MAPT
UniProt ID	TAU_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--